



Counting chambers Zwerver



Plankton
Zwerver

TECHNICAL DATA

Material	Plexi glass, glass, stainless steel
Dimensions of emptying chassis	56 x 120 mm
Inner diameter of sedimentation cylinder	24 mm
Dimensions of the chamber	plastic: 50 x 56 metal ring diameter = 44 mm
Height of sedimentation cylinders	10 ml = 17 mm 25 ml = 50 mm
Bottom glass	thickness = 0,17 mm diameter = 35 mm



Cleaning:
Normally just rinse
with water.

From time to time a
gentle wash with a
soft brush and
detergent.

**DON'T USE
HOT WATER
OR
ALCOHOL.**

How to use?
See videos for
[filling](#)
[emptying](#)